



2X Pfu mix

Cat No	Pack size	conc
PFM0500	5 ml	2x
PFM2500	5x5ml	2x

Description :*Pfu* DNA polymerase is a thermostable enzyme isolated from *Pyrococcus furiosus*. The enzyme replicates DNA at 75°C, catalyzing the polymerization of nucleotides into duplex DNA in the 5'-3' direction. *Pfu* DNA polymerase possesses 3'-5' exonuclease (proofreading) activity. Base misinsertions that may occur during polymerization are rapidly excised by the proofreading activity of the polymerase. *Pfu* DNA polymerase is recommended for use in PCR and primer extension reactions that require high-fidelity synthesis. *Pfu* DNA polymerase-generated PCR fragments are blunt-ended.

2X Pfu mix is optimized mixture contain of Pfu polymerase, reaction buffer, dNTP And enhancer as 2-fold concentration. 2x Pfu mix is designed to allow the user for quick ,easy preparation of reaction mixture. The 2x Pfu mix can be amplification PCR products up to 3-5 kb.

storage conditions: long time at -20°C
short time at 4 °C

Template

2 x Pfu mix is suitable for amplifying targets up to 3 kb from the following templates:

Genomic DNA: 10–200 ng

Plasmid DNA : 1–5 ng

cDNA : ~100 ng starting total RNA

Primers

Use 0.3 μ M per primer as a general starting point. For larger amounts of template (e.g., 200 ng genomic DNA), increasing the concentration up to 0.5 μ M per primer may improve yield.

Annealing Temperature

The annealing temperature is slightly higher than with typical PCR. The optimal annealing temperature should be ~2°C lower than the T_m of the primers used. A range of 58–68°C is recommended.

Extension Time: As little as 60 seconds per kb is suitable for most targets. Use up to 120 seconds per kb for maximum yield.

PCR Protocol:

1. Thaw the 2x Pfu mix at room temperature. Vortex the 2x Pfu mix and then spin it briefly in a micro centrifuge to collect the material in the bottom of the tube.
2. Prepare one of the following reaction mixes on ice:

Component	Volume
2x Pfu mix	12.5 ul
Primer1 (20 pmol)	1-2 ul
Primer2 (20 pmol)	1-2 ul
template	1-10 ul
ddH ₂ O	Up to 25 ul
Total	25 ul

3. If necessary you can scale up your volume

1. Program the thermal cycler as follows:

Step	Temperature	Time	Cycle
Initial denaturation	94-96°C	0.5-2mins	1
Denaturation	94-96°C	0.2-2mins	15-30
Annealing	50-68	0.2-2mins	
Extension	68-75	2min/1kb	
Final extension	68-75	1-10mins	1

Step

After cycling, maintain the reaction at 4°C. Samples can be stored at -20°C until use.

Analyze products using standard agarose gel electrophoresis.

ZymesetBiotek

Address: 4F,No.72-3 Lane 296 Sec 6 Min-Chung East Rd. 114 Taipei (Taiwan)

phone☎886)2 26347345 Fax: (886)2 26347346 email: max@zymeset.com