



## 2x PCR mix

Cat No	Pack size	conc
PM0500	5 ml	2x
PM2500	5x5ml	2x

PCR mixes (2x) contain all components for PCR, only template DNA and primers should be added, volume adjusted by water if necessary and PCR can be started.

**Description:** 2x PCR mix is optimized mixture contain of Taq enzyme, reaction buffer, dNTP and enhancer as 2-fold concentration. 2x PCR mix is designed to allow the user for quick and easy preparation of reaction mixture. The PCR mix may be amplification PCR products to 3 kb and the products can be directly cloning into T-vector.

**storage conditions:** long time at -20°C  
short time at 4 °C

### Template

2 x PCR mix is suitable for amplifying targets up to 3 kb from the following templates:

Genomic DNA: 10–200 ng

Plasmid DNA : 1–5 ng

cDNA : ~100 ng starting total RNA

### Primers

Use 0.3  $\mu$ M per primer as a general starting point. For larger amounts of template (e.g., 200 ng genomic DNA), increasing the concentration up to 0.5  $\mu$ M per primer may improve yield.

### Annealing Temperature

The annealing temperature is slightly higher than with typical PCR. The optimal annealing temperature should be ~2°C lower than the T<sub>m</sub> of the primers used. A range of 58–68°C is recommended.

**Extension Time:** As little as 30 seconds per kb is suitable for most targets. Use up to 60 seconds per kb for maximum yield.

### PCR Protocol:

1. thaw the 2x PCR mix at room temperature. Vortex the 2x PCR mix and then spin it briefly in a micro centrifuge to collect the material in the bottom of the tube.
2. Prepare one of the following reaction mixes on ice:

Component	Volume
2x PCR mix	12.5 ul
Primer1 (20 pmol)	1-2 ul
Primer2 (20 pmol)	1-2 ul
template	1-10 ul
ddH <sub>2</sub> O	Up to 25 ul
Total	25 ul

3. If necessary you can scale up your volume

### 1. Program the thermal cycler as follows:

Step	Temperature	Time	Cycle
Initial denaturation	94-96°C	0.5-2mins	1
Denaturation	94-96°C	0.2-2mins	15-30
Annealing	50-68	0.2-2mins	
Extension	68-75	1min/1kb	
Final extension	68-75	1-10mins	1

### Step

After cycling, maintain the reaction at 4°C. Samples can be stored at -20°C until use.

Analyze products using standard agarose gel electrophoresis.