



## T4 DNA ligase

Cat No	Pack size	conc
TL0500	500 U	5U/ul
TL2500	2500 U	5U/ul

### Description:

T4 DNA Ligase catalyzes the formation of a phosphodiester bonds between 5'-phosphate and 3'-hydroxyl ends in double strand DNA or RNA with blunt or cohesive-end termini.

### Applications

Joining double-stranded DNA with cohesive or blunt termini and 3'T-A overhangs.

**Storage: -20°C**

### Storage Buffer

20mM Tris-HCl (pH 7.5), 1mM DTT, 50mM KCl, 0.1mM EDTA and 50% glycerol

### 10X Ligation Buffer

Contain with Tris-HCl, MgCl<sub>2</sub>, DTT, ATP and PEG(pH 7.8 at 25oC)

**Inactivation:** By heating at 65C for 10min.

### Unit definition

0.01 Weiss unit is defined as the amount of enzyme required to catalyze the ligation of greater than 95% of the Hind III fragments of 1mg of Lambda DNA at 16°C in 20 minutes.

**Purity:** No detectable DNA and nuclease activities.

### Additional information

If the reaction is blunt end or T-A ligation must have more units of T4 DNA ligase. T4 DNA ligase is strongly inhibited by NaCl or KCl if the concentration exceeds 200mM.

### Protocol

We recommend the molar ratio of vector/insert is 3/1 to 1/3. For example: cohesive end ligation (vector is 3 kb , insert is 600bp)

Ratio of vector/insert	1/1	1/3	3/1
vector DNA	100ng	100ng	100ng
insert DNA	20ng	60ng	~7ng
10Xligation buffer	1ul	1ul	1ul
T4 DNA ligase	1ul	1ul	1ul
add H <sub>2</sub> O to final volume	10 ul		
incubate the reaction at RT for 5-30 mins			